

Lab::Measurement

**Instrumentation control with Perl –
The Next Generation**


Lab::Measurement is written in Perl

- Written in Perl (and some low-level stuff in C)
- To be used from programs written in Perl

- Perl: interpreted scripting language
- Runs on almost every OS
- Extremely good in reading data files, manipulating data, etc.
- Allows to write quick and dirty scripts that get the job done
- Also allows to write clean fullsize programs

=> ideal for experimental physics

Modular architecture of the core stack



... looks complicated, but:

Using Lab::Measurement

```
#!/usr/bin/perl

use strict;
use Lab::Instrument::HP34401A;

my $hp_gpib=$ARGV[0];

print "Reading HP34401A at GPIB address $hp_gpib\n";

my $hp=new Lab::Instrument::HP34401A(
    connection_type=>'LinuxGPIB',
    gpib_address => $hp_gpib,
    gpib_board=>0,
);

my $volt=$hp->$get_voltage_dc(10,0.00001);

print "Result: $volt V\n";
```

A real measurement!

- How about doing a gate sweep, 100 steps from 0V to 1V?

```
#!/usr/bin/perl

use Lab::Instrument::HP34401A;
use Lab::Instrument::Yokogawa7651;

# Connect to gate voltage source
my $yoko=new Lab::Instrument::Yokogawa7651({
    connection_type => 'LinuxGPIB',
    gpib_address    => 14,
    gate_protect    => 0,
});
# Connect to multimeter
my $hp=new Lab::Instrument::HP34401A(0, 21);

for (my $volt=0; $volt<=1; $volt=$volt+0.01) {
    # set Yokogawa
    $yoko->set_voltage($volt);

    # wait a second
    sleep(1);

    # read multimeter
    my $vmeas=$hp->read_voltage_dc();

    # print values
    print $volt,"\t",$vmeas,"\n";
}
```


Complex measurements possible

- Not just one input, multi-dimensional measurements
- Read out as many instruments as you like, as often as you like
- Combine with monitoring of setup (pressures, temperatures, ...)
- Low-level aspects (bus, connection, interface)
 - Want to access as many devices as possible
 - Linux / Windows
 - Different device drivers, different hardware interfaces
 - Thread-safe
- High-level aspects (your script, measurement)
 - Want to make measuring as comfortable as possible
 - Access different hardware with same / similar commands
 - Keep track of metadata
 - Plotting during measurement, easy re-plotting afterwards

High-level classes: Lab::Measurement & friends

- Provide additional tools to write better measurement scripts
- Store **metadata** alongside **data**
 - date and time
 - settings of additional instruments
 - ratio of voltage divider
 - color of the shirt you are wearing
 - everything that might be important for a later interpretation of the data
- Don't repeat yourself
 - Use above collected information automatically
 - Automatically plot data with correct axes, scaling, labels etc.

High-level classes: metadata philosophy



Axes

- Axis "bias voltage" (C1)
- Axis "diode current" (C2 / R)
- Axis "diode resistance" ($R * (C1 / C2 - 1)$)
- unit
- expression
- description...

Constants

- $R = 1000 \Omega$

Columns

- C1: V_{bias}
- C2: V_{resistor}
- unit
- description...

1.5	+8.19300500E-01
1.4	+7.23413000E-01
1.3	+6.28083900E-01
1.5	+8.19300500E-01

Plots

- Plot "diode current" axis "diode current" over axis "bias voltage"
- Plot "diode resistance" axis "diode resistance" over axis "bias voltage"
- logscale
- grid
- ranges...

Using Lab::Measurement

```
my $measurement=new Lab::Measurement(  
  sample      => $sample,  
  title       => $title,  
  filename_base => 'zener_kennlinie',  
  description  => $comment,  
  
  live_plot   => 'diode current',  
  
  constants  => [  
    {  
      'name'      => 'R',  
      'value'     => '1000',  
    },  
  ],  
  columns    => [  
    {  
      'unit'      => 'V',  
      'label'     => 'V_{bias}',  
      'description' => 'Bias Voltage',  
    },  
    {  
      'unit'      => 'V',  
      'label'     => 'Amplifier output',  
      'description' => 'Voltage drop on serial resistor',  
    }  
  ],  
  axes       => [  
    {  
      'unit'      => 'V',  
      'expression' => '$C0',  
      'label'     => 'V_{bias}',  
      'description' => 'Bias voltage',  
      'min'       => ($start_voltage < $end_voltage)  
        ? $start_voltage  
        : $end_voltage,  
      'max'       => ($start_voltage < $end_voltage)
```

```
$measurement->start_block();  
  
for (  
  my $volt = $start_voltage;  
  ($volt-$end_voltage)/$step < 0.5;  
  $volt += $step  
) {  
  $knick->set_voltage($volt);  
  sleep(1);  
  my $meas = $hp->get_voltage_dc(10,0.0001);  
  $measurement->log_line($volt,$meas);  
}  
  
my $meta = $measurement->finish_measurement();
```

Run measurement!

Describe measurement!

Result: two files per measurement

„.dat“: your measurement data, in simple text format

```
# comment
# comment
x1      x2      x3      y1      y2      y3
x1      x2      x3      y1      y2      y3
...
```

„.meta“: all the additional metadata in computer-readable xml

- can be read out and evaluated again
- can be used to automatically replot the data for example


Other high-level features

- **“Gate protect”** safety mechanism
 - Makes sure that no voltage is changed too fast
 - Big voltage steps are automatically split into small, slow steps
- **Date/time** handling
 - Date/time column can contain timestamp for every data point
 - Plot data as function of time
- Measurements with **higher dimensionality**
 - Each trace/sweep/line is a “block”
 - Two-dimensional plots, selections of traces, etc.

Utilities: `plotter.pl`, `make_overview.pl`

- **`plotter.pl`:**
 - Reads measurement file, list available plots (axes etc)
 - Creates postscript or pdf output
- **`make_overview.pl`:**
 - Reads all measurements e.g. in a directory
 - Generates a postscript or pdf file with plots for each measurement, including the metadata (i.e. constants, parameters, color of your shirt) (via LaTeX)
 - Great for completing your lab book

Modular architecture of the core stack



Internal architecture

- Divided into several layers
- **Lab::Measurement** is the highest abstraction layer. Provides support for writing good measurement scripts. Offers means of saving data and meta information to disk, plotting data, etc.
- **Lab::Instrument** package makes communication with instruments easier by silently handling the protocol involved
- **Lab::Connection** and **Lab::Bus** handle communication with the hardware and encapsulate the actual device drivers.
- The lowermost layer is given by the hardware driver and its Perl binding. Several backends are supported, e.g. NI-VISA via **Lab::VISA** or **LinuxGPIB**.

VISA, GPIB, etc.

- Instruments can be connected in various ways: Serial port, GPIB, VXI, TCP/IP, USBtm, ...
- **GPIB** (hardware and software)
 - GPIB (IEEE488): Standard by Hewlett-Packard
 - Physical layer IEEE488.1
 - Command layer IEEE488.2
 - SCPI (Standard Commands for Programmable Instruments)
- **VISA** (software)
 - Virtual Instrument Software Architecture
 - VXI, GPIB, serial, or computer-based instruments
 - NI-VISA library is one implementation of the VISA standard

Open source, free software

- Open source, free software
- License: same as Perl (Artistic / GPL-2)
- Homepage: <http://www.labmeasurement.de/>
- Releases available from CPAN
- Development code and history on Gitorious
<https://www.gitorious.org/lab-measurement/lab>
- Contributors welcome!

